



UNLIMITED MOBILE PERFORMANCE TESTING

Designed for outdoor performance diagnostics, the MetaMax 3B allows maximum wireless ranges of more than 1,000 m making spiroergometry extremely independent and flexible. The MetaMax 3B can be used in nearly any kind of sport and is designed for the use under extreme environment conditions.

Test athletes and professionals under real-time training and competitive conditions in order to obtain much more significant test results when comparing data collected under labor intensive environments.



METAMAX 3B PACKAGE

- Standard Feature
- Optional Feature




		METAMAX 3B PACKAGE
SOFTWARE	MetaSoft Studio Basic License	●
	Performance Assessment	●
	Exercise Scheduler	●
	MetaControl 3000 Interface	○
	ECG Interface	○
	3-channel ECG Interface	○
	Spirometry	○
	Metabolism for CPET	●
	Basal Resting Metabolic Rate	○
	Pre-operative Risk Assessment	●
	HL7 Interface	○
Network License Extension	○	
SPECIFICATIONS	Type	Bluetooth
	Frequency	2.4 GHz
	Heart Rate	Polar H7
	Range	1,000 m
	Interface Telemetry	USB 2.0
	PC Interface	Bluetooth 2.1
	Unit Weight	1.3 lb. (0.59 kg) (without battery)
	Width	4.3" (10.9 cm)
	Length	4.7" (11.9 cm)
	Height	1.8" (4.6 cm)
	Power Supply	Battery (6 Hour Max) (Back up included)
	Warm-Up Time	approx. 30 min
	Temperature	14° to max. +104 F
	Ambient Pressure	500 - 1,050 mbar
	Relative Humidity	0 - 99% (non-condensing)
	Measuring Method	Breath-by-breath
	Sample Line .60m/2m/3m	.60 m ● / 2 m ○ / 3 m ○
	Face Mask (S, M, L)	2 S / 2 M / 2 L
	MaxSport Jacket (S, L)	●
	Remote Control	●
	Long Range Bluetooth Receiver	●
	GPS Kit	●
	Fast Charger	○
	Back Carrying System	●
	Calibration Kit	●
	Stationary Kit	○
Warranty	1 Year Non-Wear Items 1 Year Wear Items	



TOTAL SPIROERGOMETRY SOLUTION

- Maximum wireless ranges of more than 1,000 m
- GPS measuring for the identification of velocity and altitude change
- Applicable for all kind of sport for optional back or breast wearing system
- Touch screen remote allows you to see and carry out essential data during a field test without PC or laptop



ANYWHERE & EVERYWHERE

- Designed for the use under extreme environment conditions
- Test athletes under real training and competition conditions
- Using the optional stationary kit, the MetaMax 3B can be enhanced to a full stationary system and is 100% compatible to spirometry test setup.
- Stable connections makes connecting cables easy

*Specifications subject to change without written notice.
Warranties outside the U.S. may vary - Please contact your dealer for details.*

**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT
VISIT WWW.WOODWAY.COM OR CALL 800.WOODWAY TO TALK TO A REPRESENTATIVE.**