

THE POSTURAL LINE

The Balance Trunk is a dynamic trunk balance training and assessment system performed in the sitting position. Paired with the very simple and intuitive TecnoBody software it is easy for fitness and rehabilitation professionals to evaluate and assess each client with precise and detailed data to use to measure progress and create powerful and specific training plans and programs.




BALANCE TRUNK E

- Standard Feature
- Optional Feature




BALANCE TRUNK E		
USER INTERFACE	16" Touchscreen Display	●
	Adjustable Display Height	●
	Processor Frequency	1.8 GHz
	RAM Memory	4GB
	Wireless Keyboard & Mouse	●
	Wireless Printer	○
	TecnoBody Management System	○
SPECIFICATIONS	User Weight Capacity	330 lb.
	Dynamic Platform	
	Support Platform	
	Unit Weight	176.4 lb. (shipping weight 225 lb.)
	Width	31.9" (81 cm)
	Length	63" (160 cm)
	Height	51.1" (130 cm)
	Trunk Sensor	○
	Wrapping Upper Supports	
	Safety Step	
	Trunk Balance Seat	●
	110/240 Vac 20 Amp	●
	Warranty	1 Year All Components 1 Year Labor Lifetime Software

*Specifications subject to change without written notice.
Warranties outside the U.S. may vary - Please contact your dealer for details.*



SOFTWARE

- Assessment modules including the star test & proprioceptive
- Training modules including dynamic, skill and equilibrium
- Patient data file
- Games including equilibrium and flight simulator



FUNCTIONALITY

- Dynamic seat balance
- Sitting position with adjustable seat height
- 10 different selectable levels of instability
- Seat tilt measurement: +/- 15° (maximum resolution of 0.1° and 0.5° precision)
- Can be repositioned with a wireless connection to a PC
- Ability to connect to TecnoBody Management System by TecnoBody key



**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT
VISIT WWW.WOODWAY.COM OR CALL 800.WOODWAY TO TALK TO A REPRESENTATIVE.**