



## THE FUNCTIONAL LINE

ISO-FREE is the world's first functional training system for free exercises incorporating the ability to recognize movements and balance for controlled and effective training. Engineered for sports medicine, rehabilitation and occupational therapy the ISO-FREE can accommodate users on a motion analysis and force platform to record, check and/or correct movements, balance and posture.



# ISO-FREE


- Standard Feature
- Optional Feature



USER INTERFACE

SPECIFICATIONS


	ISO-FREE
23" Touchscreen Display	●
27" Touchscreen Display	○
Processor Frequency	2.5 GHz
RAM Memory	8GB
Wireless Keyboard & Mouse	●
Wireless Printer	○
TecnoBody Management System	●
User Weight Capacity	660 lb.
Operative Surface Dimensions	39" x 39" (100 cm x 100 cm)
Access Height	5.9" (15 cm)
Unit Weight	375 lb. (shipping weight 450 lb.)
Width	51" (130 cm)
Length	96" (244 cm)
Height	75" (191 cm)
Precision	7 oz
110 Vac 20 Amp	●
Warranty	1 Year All Components 1 Year Labor Lifetime Software



## SOFTWARE

- Multi-system approach to fitness and rehabilitation
- Capabilities for assessing and managing controlled parameters simultaneously with real-time feedback.
- Multiple of training and assessment modules as well as virtual games for functional training of the upper and lower limbs

*Specifications subject to change without written notice.  
Warranties outside the U.S. may vary - Please contact your dealer for details.*



## FEATURES

- Test modules with balance and posture control as well as jump analysis, reaction speed and other tests
- Training modules for isotonic workouts
- Create powerful and specific training plans and programs with an automatic and autonomous execution
- Full integration with TecnoBody Management System and patient data files
- 23" PC Control touchscreen for real-time feedback
- 3D camera for motion capture
- Virtual games for functional training for upper and lower limbs

**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT  
VISIT [WWW.WOODWAY.COM](http://WWW.WOODWAY.COM) OR CALL 800.WOODWAY TO TALK TO A REPRESENTATIVE.**